

**To use this Pulse Oximeter and Heart Rate Monitor**, just insert your finger like the diagram indicated below\*, then press the power button and the Pulse Oximeter and Heart Rate Monitor will display your heart beat p/m, the percentage of oxygen in your blood and also displays an LED indicator bar that increases and decreases by syncing with the rhythm of your heart beat. The bar graph serves the purpose of letting you know if you've got a weak or strong hear based on the level of the graph. You get all this without needing to go to the hospital, see a doctor or have any invasive procedure.

This Pulse Oximeter and Heart Rate Monitor is suitable for physicians, home doctors (mom), medical students, healthcare workers, athletes, pilots, community sports centers etc. Demand for a lightweight, inexpensive monitor for measuring SPO2 oxygen saturation and heart rate during physically active sports and high-altitude activities is in huge demand, so we present this **Pulse Oximeter and Heart Rate Monitor** at a wholesale price and directly available to B2C and B2B consumers.

#### **At a Glance...**

- Oximeter
- Measure SPo2 level
- Measure Heart rate p/m
- LED display, simple and effective
- Excellent for high altitudes sports